



1 | Half roll



2 | Roll



3 | Bobath roll



✓ Choose from 2 cover colours

## Sport-Thieme support cushions

Very good postural characteristics due to the use of high-quality polyether foam. Hygienic, easily cleaned

and disinfected. Covered with skin-friendly, synthetic leather (100% polychloride), with zip.

Version	White	Blue
<b>1   Half roll</b>		
40x12x6 cm	71 219 6103	71 219 6116
50x18x9 cm	71 219 6204	71 219 6217
50x24x12 cm	71 262 6909	71 262 6912
60x25x12.5 cm	71 219 6305	71 219 6318
<b>2   Roll</b>		
40x12 cm	71 219 6608	71 219 6611
50x18 cm	71 219 6709	71 219 6712
50x24 cm	71 262 7003	71 262 7016
60x15 cm	71 262 5401	71 262 5414
60x25 cm	71 219 6800	71 219 6813
<b>3   Bobath roll</b>		
100x20 cm	71 219 5706	71 219 5719
100x30 cm	71 219 5807	71 219 5810
100x40 cm	-	71 219 5865
100x50 cm	-	71 219 5878
120x40 cm	71 219 5908	71 219 5911
150x50 cm	71 219 6002	71 219 6015
<b>4   Head/tummy support (curved)</b>		
40x30x9 cm	71 219 6406	71 219 6419
<b>5   Cuboid/posture support</b>		
25x15x5 cm	71 262 5502	71 262 5515
50x25x5 cm	71 262 7104	71 262 7117
50x40x20 cm	71 219 6901	71 219 6914
50x45x40 cm	71 219 7005	71 219 7018
60x50x40 cm	71 219 7106	71 219 7119
<b>6   Wedge</b>		
45x35x12/2 cm	71 219 6507	71 219 6510
<b>7   Leg support</b>		
45x42x20 cm	71 219 7207	71 219 7210

4 | Pillow



5 | Cuboid



6 | Wedge



7 | Leg support



## 8 | Majimo Cushion

Thanks to its perfect orthopaedic and ergonomic shape, the cushion is excellent as a positional aid. Using the 3-sided cushion, patients can be easily and safely re-positioned in a way that saves time and effort. The special shape of the cushion makes it easy to immediately change a patient's position simply by pulling on one end of the cushion. The patient no longer needs to be lifted for every change in position. The cushion retains its shape and has no seams on the long sides that face the body: this also prevents the material from being displaced. In combination with the ability to change positions quickly and more often, this prevents the development of pressure sores. The cushion is suitable for people with allergies and is particularly hygienic. Ticking: fibre composite made of 50% TENCEL and 50% polyester, filling: hollow polyester fibre not treated with silicone (Dacron 95). Washable at up to 95°C, can be tumble dried. Available in 3 sizes for different types of positioning. 'Small', side length: 52 cm, 'Medium', side length: 72 cm, 'Large', side length: 92 cm.

71 270 8805 Small Each  
71 270 8818 Medium Each  
71 270 8821 Large Each

## 9 | Majimo Cushion Covers


In 3 different sizes suitable for use with the Majimo cushion. Available versions: fabric cover or waterproof protective cover. The fabric cover is made of a blended fibre consisting of 50% TENCEL and 50% polyester. The waterproof cover is coated with laminated membrane made of polyurethane that prevents liquids, bacteria, viruses and mites from penetrating the cushion.

### Fabric cover

71 270 8906 Small Set of 2  
71 270 8919 Medium Set of 2  
71 270 8922 Large Set of 2

### Waterproof cover

71 270 8935 Small Each  
71 270 8948 Medium Each  
71 270 8951 Large Each



**New!**

## 10 | New! Theraline 'Pearlfusion' Neck Cushion and Pillow

The perfect support for the shoulders and neck: the neck cushion and pillow is ideal for those that sleep on their sides or backs and provides optimal comfort. The upper layer of the cushion is made of viscoelastic foam and the base layer is made of micro pearls. This allows the cushion to adapt and perfectly fit your head and neck, providing excellent stability. Outer cover: 100% cotton, inner cover: 87% cotton, 13% elastane. Filling: 100% polystyrene micro pearls (80 g/l). Filling: foam core (polyurethane). Foam core cannot be washed. Covers can be washed at up to 60°C – can be tumble dried.

• Standard width: 50x32x12 cm, 2 kg  
• Comfort width: 70x32x12 cm, 3.8 kg  
71 286 4600 Standard width Each  
71 286 4613 Comfort width Each